



General Instructions for Culturing Bovine Renal Arterial Endothelial Cells (BRAEC)

Open all the packages immediately upon arrival and examine each component for shipping damage. Notify Cell Applications, Inc. or your distributor immediately if there is any problem.

I. STORAGE

A. CRYOPRESERVED VIALS (B900-05)

Store the cryovials in a liquid nitrogen storage tank immediately upon arrival.

B. PROLIFERATING FLASKS (B901-25, -75)

1. Examine under a microscope to check if all the cells are attached to the bottom of the flask. If not, notify CAI or your distributor immediately.
2. Decontaminate the exterior of the flask with 70% alcohol.
3. Place the sealed flask in a 37°C, 5% CO₂ humidified incubator for 2 hours as shipped.
4. In a sterile Biological Safety Cabinet, open the cap of the flask very slowly and carefully.
5. Remove the Transport Medium by aspiration. Add fresh Growth Medium: 5 ml for a T-25 flask and 15 ml for a T-75 flask.
6. Place the flask in a 37°C, 5% CO₂ humidified incubator with loosened cap to allow gas exchange.
7. Change medium every other day.

C. GROWTH MEDIUM (B211-500)

Store the Growth Medium at 4°C in the dark immediately upon arrival.

D. REAGENT KIT (090K)

Store at -20°C immediately upon arrival.
Store at 4°C after thawing.

II. PREPARATION FOR CULTURING

1. Make sure the Class II Biological Safety Cabinet, with HEPA filtered laminar airflow, is in proper working condition.
2. Clean the Biological Safety Cabinet with 70% alcohol to ensure it is sterile.
3. Turn the Biological Safety Cabinet blower on for 10 min. before cell culture work.
4. Make sure all serological pipettes, pipette tips and reagent solutions are sterile.
5. Follow the standard sterilization technique and safety rules:
 1. Do not pipette with mouth.
 2. Always wear gloves and safety glasses when working with animal cells even though all the animals have been inspected by USDA.
 3. Handle all cell culture work in a sterile hood.

III. CULTURING BRAEC

1. PREPARING CELL CULTURE FLASKS FOR CULTURING BRAEC

1. Take the Attachment Factor Solution and Bovine Endothelial Cell Growth Medium from the refrigerator. Decontaminate the bottles with 70% alcohol in a sterile hood.
2. Swirl the Attachment Factor Solution bottle a few times to form a homogenous solution. If the solution is gelled, warm up in a 56°C water bath for 10 minutes and swirl the bottle to mix the solution.
3. Add 7.5 ml of Attachment Factor Solution* to a T-75 flask and rock gently to distribute the solution evenly to cover the whole culture surface.
4. Coat the culture ware for 30 minutes at 37°C or 2 hours

European Collection of Cell Cultures
Health Protection Agency
Centre for Emergency Preparedness and Response
Porton Down
Salisbury
SP4 0JG, UK
Tel : +44 (0) 1980 612512
Fax: +44 (0) 1980 611315
Email: ecacc@hpa.org.uk
Website: www.ecacc.org.uk

Cell Applications Inc.
5820 Oberlin Drive
Suite 101
San Diego
CA 92121, USA

Tel: (858) 453-0848 Toll-Free: (800) 645-0848
Fax: (858) 453-2862
Email: info@cellapplications.com
Web: www.cellapplications.com

(overnight is OK) at room temperature.

5. Remove the Attachment Factor Solution by aspiration in a sterile hood.
6. The coated flask can be used immediately or stored at 4°C for up to one month.
7. Prepare the coated flask for subculturing by pipetting 12 ml of Bovine Endothelial Growth Medium into this coated T-75 flask and wait for seeding.

*The coating concentration is 1 ml per 10 cm² surface area of culture ware.

-2.5 ml for coating a T-25 flask or a 60mm culture dish.

-7.5 ml for coating a T-75 flask or a 100mm culture dish.

B. THAWING AND PLATING BRAEC

1. Remove the cryopreserved BRAEC vial from the liquid nitrogen storage tank using proper protection for your eyes and hands.
2. Turn the vial cap a quarter turn to release any liquid nitrogen that may be trapped in the threads, then retighten the cap.
3. Thaw the cells quickly by placing the lower half of the vial in a 37°C water bath for 1-2 minutes.
4. Take the vial out of the water bath and wipe dry.
5. Decontaminate the vial exterior with 70% alcohol in a sterile Biological Safety Cabinet.
6. Remove the vial cap carefully. Do not touch the rim of the cap or the vial.
7. Resuspend the cells in the vial by gently pipetting the cells 5 times with a 2 ml pipette. Be careful not to pipette too vigorously as to cause foaming.
8. Pipette the cell suspension (1ml) from the vial into the T-75 flask that is pre-coated with Attachment Factor Solution and contains 1.2 ml of Bovine Endothelial Cell Growth Medium.
9. Cap the flask and rock gently to evenly distribute the cells.
10. Place the T-75 flask in a 37°C, 5% CO₂ humidified incubator. Loosen the cap to allow gas exchange. For best results, do not disturb the culture for 24 hours after inoculation.
11. Change the Bovine Endothelial Cell Growth Medium after 24 hours or overnight to remove all traces of DMSO.
12. Change Bovine Endothelial Cell Growth Medium every other day until the cells reach 60% confluent.
13. Double the Bovine Endothelial Cell Growth Medium volume when the culture is >60% confluent or for weekend feedings.
14. Subculture the cells when the BRAEC reach 80%

confluent.

IV. SUBCULTURING BRAEC

A. PREPARING SUBCULTURE REAGENTS

1. Remove the Subculture Reagent Kit from the -20°C freezer and thaw overnight in a refrigerator.
2. Make sure all the subculture reagents are thawed. Swirl each bottle gently several times to form homogeneous solutions.
3. Store all the subculture reagents at 4°C in the dark for future use. The activity of Trypsin/EDTA Solution will be stable for 2 weeks when stored at 4°C.
4. Aliquot Trypsin/EDTA solution and store the unused portion at -20°C if only portion of the Trypsin/EDTA is needed.

B. PREPARING CULTURE FLASK

Follow the same procedure in IIIA.

C. SUBCULTURING BRAEC

Trypsinize Cells at Room Temperature. Do Not Warm Any Reagents to 37°C.

1. Remove the medium from culture flasks by aspiration.
2. Wash the monolayer of cells with HBSS and remove the solution by aspiration.
3. Pipette 5 ml of Trypsin/EDTA Solution into the T-75 flask. Rock the flask gently to ensure the solution covers all the cells.
4. Remove 4.5 ml of the solution immediately.
5. Replace the cap tightly and monitor the trypsinization progress at room temperature under an inverted microscope. It usually takes about 1 minute for the cells to become rounded but still attached to the flask. (When rounded cells detach by itself without hitting, it means the cells are over trypsinized.)
6. Release the rounded cells from the culture surface by hitting the side of the flask against your palm until most of the cells are detached.
7. Pipette 5 ml of Trypsin Neutralizing Solution to the flask to inhibit further tryptic activity.
8. Transfer the cell suspension from the flask to a 50 ml sterile conical tube.
9. Rinse the flask with an additional 5 ml of Trypsin Neutralizing Solution and transfer the solution into the same conical tube.
10. Examine the T-75 flask under a microscope. If there are >20% cells left in the flask, repeat Steps 2-9.

European Collection of Cell Cultures
Health Protection Agency
Centre for Emergency Preparedness and Response
Porton Down
Salisbury
SP4 0JG, UK
Tel : +44 (0) 1980 612512
Fax: +44 (0) 1980 611315
Email: ecacc@hpa.org.uk
Website: www.ecacc.org.uk

Cell Applications Inc.
5820 Oberlin Drive
Suite 101
San Diego
CA 92121, USA

Tel: (858) 453-0848 Toll-Free: (800) 645-0848
Fax: (858) 453-2862
Email: info@cellapplications.com
Web: www.cellapplications.com

11. Centrifuge the conical tube at 200 x g for 5 minutes to pellet the cells.
12. Remove the conical tube. Aspirate the supernatant.
13. Flick the tip of the conical tube with your finger to loosen the pellet.
14. Resuspend the cells in 5 ml of Bovine Endothelial Cell Growth Medium by gently pipetting the cells to break up the clumps.
15. Count the cells with a hemocytometer or cell counter. Inoculate at 10,000 cells per cm² for rapid growth, or at 5,000 cells per cm² for regular subculturing.
16. Transfer the 5 ml cell suspension to a T-225 flask containing 45 ml of Bovine Endothelial Cell Growth Medium, and do 1 to 4 split thereafter.

European Collection of Cell Cultures
Health Protection Agency
Centre for Emergency Preparedness and Response
Porton Down
Salisbury
SP4 0JG, UK
Tel : +44 (0) 1980 612512
Fax: +44 (0) 1980 611315
Email: ecacc@hpa.org.uk
Website: www.ecacc.org.uk

Cell Applications Inc.
5820 Oberlin Drive
Suite 101
San Diego
CA 92121, USA

Tel: (858) 453-0848 Toll-Free: (800) 645-0848
Fax: (858) 453-2862
Email: info@cellapplications.com
Web: www.cellapplications.com